Defibrillators

- 1. Woodfold View South West Ambulance service is responsible for the defibrillator and servicing
- 2. Telephone Box Corscombe Friends of the First Responders are responsible for, along with the Parish Council.
- 3. Village Shop machine South West Ambulance service is responsible for the defibrillator and servicing
- 4. Halstock Village Hall machine Friends of the First Responders are responsible for, along with the Parish Council.
- 5. West Chelborough Friends of the First Responders are responsible for, along with the Parish Council.

Defibrillator Report

5th June 2023

- 1. Inspection at around 9 pm
 - a. Telephone Box Corscombe fan off.
 - i. Machine ready for action
 - ii. Dick Thorne (Chris C as deputy) looks after this machine.
 - b. Woodfold View fan off but defib (or cabinet) beeping.
 - i. Battery in defib flat
 - ii. Reported to SWASFT by Geoff Williams but no response as yet.
 - iii. Geoff Williams looks after this Machine
 - c. Village Hall Halstock fan on
 - Returned at 10 pm fan off.
 - ii. Machine ready for action
 - iii. Peter Brinck looks after this machine.
 - iv. Suggested that he recruits a deputy
 - d. Shop Halstock no fan fitted to cabinet
 - i. Machine ready for action
 - ii. Old style cabinet which is under cover so shaded

- iii. Gill Clark looks after this unit
- e. West Chelborough fan off
 - i. Machine ready for action
 - ii. Susan Soames looks after this machine.
 - iii. Suggested Susan recruits a deputy.
- 2. Design fault with new cabinets

 - a. Fan air inlet / outlet partially blocked by defib.i. Make a small stand for defib so air vent is clear.
 - ii. Have a reflector vest for slipping over the defib case.

Many thanks to Chris C for this update